

# GreenCare<sup>®</sup>

MEASURED MOBILITY



## GENERAL PURPOSE CHAIR

# BRITISH BUILT WHEELCHAIRS

Using the excellent components from the Greencare bespoke wheelchair range this general purpose folding wheelchair combines strength and comfort with light weight. It comes in three standard seat widths as either an attendant propelled or for the more able, self-propelling wheels.

The Greencare general purpose chair offers additional configurability using the range of Greencare accessories.

## STANDARD FEATURES

- Seat sizes 15" wide and 15" deep, 17" wide and 17" deep and 18" wide and 17" deep all on single or double braced frames. Suitable for occupants to 125KG.
- Seat to ground a standard height of 17½".
- Rear wheel of 24" (610mm) for self-propelling or 16" (405mm) for carer assisted chair (Pic 1).
- Armrest height of 10" (250mm) complete with arm pad.
- Backrest height of 19" (480mm).
- Vehicle tie down points for transportation (Pic 2).
- Lap belt fitted as standard (Pic 3).
- Cushioned non-sag fire retardant fabric for seat and back rest ensures comfort.
- Large front castors make pushing a little easier.
- All tyres made by Greentyre and are long life, non marking and puncture proof.
- Scratch resistant easy clean anodised non painted frame.
- Greencare wheelchairs are crash tested to ISO 7176-19.



## IT'S THE LITTLE THINGS THAT MAKE GREENCARE SPECIAL

Chairs are built by Greencare in the United Kingdom and can be customised to meet individual needs. Using quality components designed by the finest British engineers Greencare chairs can be matched by no other as they are the most versatile and uniquely configurable chairs around.

Important dimensions such as seat depths, seat to ground height other than those shown as standard are possible.

You won't be alone when ordering a wheelchair from Greencare as a trained Customer Service Team are on standby to help answer any questions you may have.



Greencare wheelchairs are environmentally friendly, they use the finest aircraft grade aluminium, have a natural anodised finish that minimises unsightly scuff marks commonly experienced on painted surfaces.

We take pride in making long lasting comfortable wheelchairs that are good for the occupant but also work in the environment in which they have to operate, taking into account tight spaces, ramps, small doors and difficult surfaces.

Greencare - a world class product.



# ACCESSORIES

Standard ●  
Optional ○


## SEAT



**Lapstrap**  
With airline type buckle.




**Thoracic support**  
Adjusts for height, width, length on non folding backs.



**Extended handles**  
To aid fitting hub brake levers or other controls.



**Extended seat**  
Brings occupant forward, requires outrigger footrests.



**Personalise**  
With coloured tyres, signwriting or green side panels.



**Outrigged frame**  
Outrigged frame that increases front stability.


## WHEELS



**Spoke guards**  
In clear or black plastic.



**Greenrim ergonomic design for easy push**  
Warm to touch and scuff free - 24" only.



**Brake lever extension**  
Fold down lever to aid side transfer.



**One arm drive**  
For 24" wheels please state drive side.



**Super heavy duty wheel**  
24" and 26" wheels tested to 320KG available with hub brakes.



**Hub brakes**  
Available on all wheels sizes. Needs extended handles option.

## FRAME




**Reinforced back post**  
Best used when backpost in forward position.



**Crutch holder**  
Left or right fixing.



**Oxygen carrier**  
Uses rear backpost receiver for fixing.



**IV pole**  
1.3m above seat rails. State left or right mount.



**Drop frame**  
Drop frame for ultra-low seat to ground.



**Reinforced side panels**  
Aluminium rigid construction.

## ARMREST



**Complete tray**  
With soft edges and cushioned front.




**Therapy tray**  
Non 'handed' cushioned tray with positioning options.



**Universal holder**  
For keeping close handy items, water, TV remote etc.



**Tray mount bracket**  
Unique self centering brackets taking the fuss out of tray mounting.



**Extended armrest pad**  
Extra long arm pad.



**Height adjustable armrests**

## LEGS



**ELR with calf pad**  
Length 23" 580mm. Can be supplied without footplate. Left or right.




**Stump supports**  
Residual limb support with cushioned pad.



**Calf straps**  
Sized to chair.



**Adductor pads**  
Adjustable pad on footrest hanger.



**High level release**  
For users who can't lean forward to release footplate catches.



**Dynamic footboard**  
Carcassone style swinging footplate.


## FEET



**Angle adjustment footplates**



**Universal footplates**  
Adjustment 80° to 115° + footplate positioning. Replaces old CAPS 5.



**Padded footplates**  
Soft cushion pad replaces foot plate.


## CUSHIONS



Viscoe elastic and high resilience foam core supplied sized to chair in 2" and 3" thickness. The fabric is washable with non slip on base.



Making the right choice of wheelchair is important, that is why we are a preferred supplier to the UK's National Health Service. We have a team on hand to help with the little decisions that will turn a great wheelchair into a perfect choice so if you are choosing a chair for you or someone you care for or love then let us help.



To order call +44 (0)1642 22 33 22  
or email [info@greencaremobility.com](mailto:info@greencaremobility.com)

## CARING AND SUPPORTING

Greencare lightweight wheelchairs are amongst the most environmentally friendly wheelchairs in the world and certainly one of the most comfortable, designed by Europe's leading wheelchair designer David Batty and manufactured in the North East England by a dedicated proud workforce.



## A VERY SPECIAL WARRANTY

Greencare chairs have a 36 month replacement part warranty on any part found to have a manufacturing defect. Technical manuals, Technical Guide's, User Guide's, assembly drawings, service information, part numbers and more are available on line from the Greencare Website [www.greencaremobility.com](http://www.greencaremobility.com) Servicing and repairs should be conducted by qualified approved repairers.

Greencare in-country support is available in Abu Dhabi, China, Dubai, Eire, Qatar and Saudi Arabia.  
© Copyright, FFT Ltd 2015 All rights reserved. Greencare Mobility, Simcox Ct, Middlesbrough TS2 1UU.  
Whilst every care was taken in preparing this leaflet E&OE.